

INDIRECT ELISA

www.ptglab.com

All steps are carried out at room temperature unless stated otherwise.
Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

-
1. **Antigen coating:**
 - a. Dilute purified antigens to a final concentration of 3.0 µg/ml in **antigen-coating buffer** and add 50 µl of diluted antigen to each well of a 96-well ELISA plate.
 - b. Carefully cover the plate with adhesive plastic and incubate at 4°C overnight.

 2. **Blocking:**
 - a. Empty the wells of **antigen-coating buffer** and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.
 - b. Add 200 µl **blocking buffer** per well to block residual protein-binding sites. Cover the plate with adhesive plastic and incubate for 1–2 h at 37°C.

 3. **Antibody incubation:**
 - a. Dilute your primary antibody of choice with **blocking buffer** in a series e.g. 1:500, 1:1000, 1:2000, 1:4000 and so on, empty the wells of **blocking buffer** and then add 100 µl of each dilution per well. Repeat in duplicate, or triplicate, for accuracy. Cover the plate with adhesive plastic and incubate for 1 h at 37°C.
 - b. Empty the wells and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.
 - c. Dilute the HRP-conjugated secondary antibody with **blocking buffer** at an optimal concentration (a dilution factor within 1:10,000-1:100,000 is recommended) and add 100 µl of secondary antibody solution to each well. Cover the plate with adhesive plastic and incubate for 1 h at 37°C.
 - d. Empty the wells and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.

 4. **Signal detection:**
 - a. Add 100 µl TMB substrate (mix equal volumes of **TMB buffer A** and **buffer B**) to each well with a multichannel pipette. Color development should peak after 15 minutes, at which time it should be stopped by adding 100 µl of 2 M H₂SO₄ per well. Read absorbance at 450 nm.

INDIRECT ELISA

www.ptglab.com

Buffers Needed

Antigen coating-buffer	For 1000 ml	PBST buffer	For 1000 ml
100 mM NaHCO ₃	8.4 g	10 mM Na ₂ HPO ₄	1.42 g
Adjust pH to 9.6		1.8 mM NaH ₂ PO ₄	0.22 g
Add ddH ₂ O to 1000 ml		140 mM NaCl	8.19 g
		0.2 % Tween 20	2 ml
		Adjust pH to 7.4	
		Add ddH ₂ O to 1000 ml	

Blocking buffer	For 100 ml
5% non-fat dry milk	5 g
Add PBST buffer to 100 ml	

TMB buffer A	For 500 ml	TMB buffer B	For 500 ml
NaAc•3H ₂ O	13.6 g	TMB (first dissolved in 3 ml DMSO)	0.15 g
Citric acid	1.6 g	EDTA-2Na	0.2 g
30% H ₂ O ₂	0.3 ml	Citric acid	0.95 g
Add ddH ₂ O to 500 ml		Glycerol	50 ml
		Add ddH ₂ O to 500 ml	