
Actin Chromobody[®]-TagRFP plasmid




For plasmid sequence, please contact info@chromotek.com

Location of features

PCMV IE: 1-589 Enhancer region: 59-465 TATA box: 554-560 Transcription start point: 583 Actin-VHH: 621-986 TagRFP: 1050-1763 Stop codon: 1761-1763 SV40 early mRNA polyadenylation signal Polyadenylation signals: 1919-1924 & 1948-1953 mRNA 3' ends: 1957 & 1969 f1 single-strand DNA origin: 2016-2471 SV40 origin of replication: 2812-2947 SV40 early promoter Enhancer (72-bp tandem repeats): 2648-2719 & 2720-2791 21-bp repeats: 2795-2815, 2816-2836 & 2838-2858 Early promoter element: 2871-2877 Kanamycin/neomycin resistance gene Neomycin phosphotransferase coding sequences: Start codon (ATG): 2999-3001; Stop codon: 3791-3793

Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal

Polyadenylation signals: 4029-4034 & 4042-4047 CoIE1 replication origin: 4325-5007



Product	Code	Size
pAC-TagRFP	acr	20 µg
Vector type	mammalian expression vector	
Reporter	TagRFP	
Reporter codon usage	mammalian	
Promoter for Chromobody®	Рсму іє	
Host cells	mammalian	
Selection	prokaryotic – kanamycin	
	eukaryotic - neomycin (G418)	
Replication	prokaryotic - pUC ori	
eukaryotic - SV40 ori		
Use	Actin Chromobody [®] -TagRFP expression in mammalian cells for non-invasive live cell visualization of endogenous	
	cytoskeleton.	

nnono

Vector description

Actin Chromobody[®]-TagRFP plasmid (pAC-TagRFP) is a mammalian expression vector encoding the cytoskeleton marker Actin-V_HH fused to red fluorescent protein TagRFP (from Evrogen). The vector allows expression Actin Chromobody[®]-TagRFP fusion protein in eukaryotic (mammalian) cells.

Chromobody[®] codon usage is optimized for high expression in mammalian cells (humanized) [Haas et al. 1996].

The vector backbone contains immediate early promoter of cytomegalovirus ($P_{CMV | E}$) for protein expression, SV40 origin for replication in mammalian cells expressing SV40 T-antigen, ColE1 origin of replication for propagation in *E. coli* and f1 origin for single-stranded DNA production. SV40 polyadenylation signals (SV40 poly A) direct proper processing of the 3'-end of the reporter mRNA.

SV40 early promoter (P_{SV40}) provides neomycin resistance gene (Neo^r) expression to select stably transfected eukaryotic cells using G418. Bacterial promoter (P) provides kanamycin resistance gene expression (Kan^r) in *E. coli.* Kan^r/Neo^r gene is linked with herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals.

Expression in mammalian cells

pAC-TagRFP vector can be transfected into mammalian cells by any known transfection method. If required, stable transformants can be selected using G418 [Gorman 1985].

Propagation in E. coli

Suitable host strains for propagation in *E. coli* include DH5alpha, HB101, XL1-Blue, and other general purpose strains. Plasmid incompatibility group is pMB1/CoIE1. The vector confers resistance to kanamycin (30 µg/ml) to *E. coli* hosts. Copy number in *E. coli* is about 500.

Note: The plasmid DNA was isolated from dam⁺-methylated *E.coli*. Therefore some restriction sites are blocked by methylation. If you wish to digest the vector using such sites you will need to transform the vector into a dam⁻ host and make fresh DNA.

Notice to Purchaser:

Chromoboly^{ma}related materials (the Products) are intended for research use only. The Products are covered by U.S. Pat. applications pending. By use of these Products, you accept the terms and conditions of the applicable End User License Agreement (EULA non-profit entities). The CMW promoter is covered under U.S. Patents 5, 168, 062 and 5, 385, 839, and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of lowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242. MATERIAL SAFETY DATA SHEET INFORMATION: To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.