


# TOLL-LIKE RECEPTOR PATHWAY

www.ptglab.com






## Related Products


| Antibody Name | Catalog Number | Type        | Applications               |
|---------------|----------------|-------------|----------------------------|
| AP1, JUN, P39 | 10024-2-AP     | Rabbit Poly | ELISA, IF, IHC, IP, WB     |
| ATF2          | 14834-3-AP     | Rabbit Poly | ELISA, IHC, WB             |
| BLK           | 10510-3-AP     | Rabbit Poly | ELISA, IHC, IP, WB         |
| CD14          | 17000-3-AP     | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| CREB1         | 12208-3-AP     | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| IKKβ/IKK1     | 15649-3-AP     | Rabbit Poly | ELISA, IF, WB              |
| IRAK1         | 10478-2-AP     | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| IRAK4         | 18221-1-AP     | Rabbit Poly | ELISA, IP, WB              |
| IRF3          | 11312-1-AP     | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| IRF7          | 22392-1-AP     | Rabbit Poly | ELISA, IHC, IP, WB         |
| IRF8          | 18977-1-AP     | Rabbit Poly | ELISA, IHC, WB             |
| JNK           | 10023-1-AP     | Rabbit Poly | ELISA, IF, IHC, IP, WB     |
| LTA           | 19157-1-AP     | Rabbit Poly | ELISA, IP, WB              |

| Antibody Name | Catalog Number | Type        | Applications               |
|---------------|----------------|-------------|----------------------------|
| MD2           | 11784-1-AP     | Rabbit Poly | ELISA, IHC, WB             |
| MKK7          | 55030-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB         |
| NFκB1         | 14220-1-AP     | Rabbit Poly | ELISA, IF, IP, WB          |
| PI3K MAPK     | 14064-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB         |
| PI3KCA/PI3K   | 20583-1-AP     | Rabbit Poly | ELISA, FC, IHC, WB         |
| RAC1          | 24072-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB         |
| RIP1          | 18443-1-AP     | Rabbit Poly | ELISA, IHC, IP, WB         |
| RIPK2         | 15366-1-AP     | Rabbit Poly | ELISA, WB                  |
| TAB1          | 14819-1-AP     | Rabbit Poly | ELISA, IP, WB              |
| TAB2          | 14410-1-AP     | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| TLR3          | 17766-1-AP     | Rabbit Poly | ELISA, IHC, IHC            |
| TLR4          | 19811-1-AP     | Rabbit Poly | ELISA, IHC, WB             |
| TLR5          | 19810-1-AP     | Rabbit Poly | ELISA, IHC, IP, WB         |

| Antibody Name | Catalog Number | Type        | Applications       |
|---------------|----------------|-------------|--------------------|
| TLR7          | 17232-1-AP     | Rabbit Poly | ELISA, IHC, WB     |
| TNKS          | 18030-1-AP     | Rabbit Poly | ELISA, IF, IP, WB  |
| TOLLIP        | 11315-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB |
| TRAF6         | 12809-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB |
| TRAM1         | 12705-1-AP     | Rabbit Poly | ELISA, IF, IHC, WB |
| UBC13/UBE2N   | 10243-1-AP     | Rabbit Poly | ELISA, IHC, IP, WB |
| UBE2V1        | 10207-2-AP     | Rabbit Poly | ELISA, WB          |

**Legend**

-  Association
-  Directly Activates
-  Inhibits

 This number shows the amount of times our antibody has been cited in a publication.