

For Research Use Only

IFITM2 Polyclonal antibody

Catalog Number: 12769-1-AP

Featured Product

39 Publications



Basic Information

Catalog Number: 12769-1-AP	GenBank Accession Number: BC009696	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop and 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 10581	Recommended Dilutions: WB 1:1000-1:4000 IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate
Source: Rabbit	Full Name: interferon induced transmembrane protein 2 (1-8D)	IHC 1:20-1:200 IF 1:20-1:200
Isotype: IgG	Calculated MW: 132 aa, 15 kDa	
Immunogen Catalog Number: AG3451	Observed MW: 15 kDa	

Applications

Tested Applications:
FC, IF, IHC, IP, WB, ELISA

Cited Applications:
FC, IF, IHC, IP, WB

Species Specificity:
human, mouse

Cited Species:
human, mouse, pig

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HepG2 cells, MCF-7 cells

IP : HepG2 cells,

IHC : human breast cancer tissue, mouse brain tissue

IF : HepG2 cells,

Background Information

IFITM2, also named as 1-8D, belongs to the CD225 family. It is an IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens, namely influenza A H1N1 virus, West Nile virus (WNV), and dengue virus, by inhibiting the early steps of replication. IFITM2 induces cell cycle arrest and mediates apoptosis by caspase activation and in p53-independent manner. It is overexpressed in colon carcinoma. IFITM2 is a novel pro-apoptotic gene that will provide new insights into the regulated cellular pathways to death. IFITM proteins are recently identified as viral restriction factors that inhibit infection mediated by the influenza A virus (IAV) hemagglutinin (HA) protein. Also they serve as important components of the innate immune system to restrict HIV-1 infection. Catalog#12769-1-AP is a rabbit polyclonal antibody produced with full-length of human IFITM2.

Notable Publications

Author	Pubmed ID	Journal	Application
Florian Wrensch	25256397	Viruses	WB
Wei Zhang	25228491	J Gen Virol	WB, IF
Jinsun Kim	33209202	Anim Cells Syst (Seoul)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 12769-1-AP (IFITM2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human breast cancer slide using 12769-1-AP (IFITM2 Antibody) at dilution of 1:50.




Immunofluorescent analysis of HepG2 cells using 12769-1-AP (IFITM2 antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG.



1X10⁶ MCF-7 cells were stained with 0.2 ug Anti-Human IFITM2 (12769-1-AP) and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:1000 (green), and 0.2 ug Control Antibody. Cells were fixed with 90% MeOH.



WB result of IFITM2 antibody (12769-1-AP; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-IFITM2 transfected HepG2 cells.



IP result of anti-IFITM2(IP:12769-1-AP, 4ug; Detection:12769-1-AP 1:10000) with HepG2 cells lysate 1360 ug.