

For Research Use Only

# ERK1/2 Polyclonal antibody

Catalog Number: 16443-1-AP **390 Publications**



## Basic Information

<b>Catalog Number:</b> 16443-1-AP	<b>GenBank Accession Number:</b> BC017832	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 400 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5594	<b>Recommended Dilutions:</b> WB 1:1000-1:4000 IHC 1:50-1:500 IF 1:10-1:100
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P28482	
<b>Isotype:</b> IgG	<b>Full Name:</b> mitogen-activated protein kinase 1	
<b>Immunogen Catalog Number:</b> AG9772	<b>Calculated MW:</b> 360 aa, 41 kDa	
	<b>Observed MW:</b> 38-44 kDa	

## Applications

<b>Tested Applications:</b> WB, IF, IHC, ELISA	<b>Positive Controls:</b> WB : NIH/3T3 cells, A431 cells, mouse kidney tissue, HEK-293 cells, HeLa cells, PC-3 cells, rat kidney tissue, mouse brain tissue, rat brain tissue IHC : human liver tissue, human kidney tissue IF : HeLa cells,
<b>Cited Applications:</b> WB, IP, IF, IHC	
<b>Species Specificity:</b> human, mouse, rat	
<b>Cited Species:</b> human, chicken, rat, sheep, mouse, pig, Eriocheir sinensis	
<b>Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0</b>	

## Background Information

ERK1 and ERK2 belongs to the protein kinase superfamily. It is involved in both the initiation and regulation of meiosis, mitosis, and postmitotic functions in differentiated cells by phosphorylating a number of transcription factors such as ELK-1. ERK1/2 catalyzed the reaction: ATP + a protein = ADP + a phosphoprotein. It is activated by tyrosine phosphorylation in response to INS and NGF. This antibody can recognize both ERK1 and ERK2 with the molecular mass of 38-44 kDa.

## Notable Publications

Author	Pubmed ID	Journal	Application
Xin-Sen Chen	36182039	Pharmacol Res	WB
YanHua Fan	36174847	Fitoterapia	WB
Weihua Xu	31702033	Mol Med Rep	WB

## Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 16443-1-AP (ERK1/2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 16443-1-AP (ERK1/2 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 16443-1-AP (ERK1/2 Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of HeLa cells using 16443-1-AP (ERK1/2 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).