

For Research Use Only

PR repeat Polyclonal antibody

Catalog Number: 23979-1-AP **24 Publications**



Basic Information

Catalog Number:

23979-1-AP

Size:

150ul , Concentration: 400 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

Dot Blot, IF, IHC, WB

Species Specificity:

human

Cited Species:

human, rat, mouse, Drosophila

Background Information

The C9orf72 "GGGGCC" repeat sequence codes five repeat peptide "GA repeat; GAGAGAGAGA", "GP repeat; GPGPGPGPG", "GR repeat; GRGRGRGRG", "AP repeat; APAPAPAPA" and "PR repeat; PRPRPRPRP". It was described previously that aggregated forms of poly-GA and poly-GP proteins do not enter the separation gel (PMID: 26374446). This antibody is used to detect the "PR repeat" sequence. This antibody detects the PRPRPRPRP peptide with dilution 1:4,500 in Elisa.

Notable Publications

Author	Pubmed ID	Journal	Application
Ian RA Mackenzie	26374446	Acta Neuropathol	WB, IHC
Junli Gao	36103513	Sci Transl Med	WB
Valeria Colicchia	36062323	FEBS Open Bio	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PR repeat Antibody (23979-1-AP) tested in ELISA. PR repeat peptide was coated onto the microplates at 1.5 µg/well and then incubated with dilution series of PR repeat antibody (23979-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.