

For Research Use Only

SARS-COV-2 NSP7 Polyclonal antibody

Catalog Number: 29263-1-AP



Basic Information

Catalog Number: 29263-1-AP	GenBank Accession Number: YP_009742614	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop;	GeneID (NCBI): 43740578	
Source: Rabbit	Full Name: ORF1ab	
Isotype: IgG		
Immunogen Catalog Number: AG30933		

Applications

Tested Applications:
ELISA

Species Specificity:
Virus

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SARS-CoV-2 Envelope protein Antibody (29263-1-AP) tested by ELISA. SARS-COV-2 NSP7 protein was coated onto microtiter plates at 1 µg/well and then incubated with a dilution series of SARS-COV-2 NSP7 protein Antibody (29263-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.