

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

P62;SQSTM1 Fusion Protein

Basic Information

Catalog Number:

Ag13131

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MASLTVKAYLLGKEDAAREIRRF5FCCSPEPEAEAEAA
AGPGPCERLLSRVAALFPAALRPGGFQAHYRDEGDGLV
AFSSDEELTMAMSYVKDDIFRIYIKEKKECRDRHPPC
AQEAPRNMVHPNVICDGCNGPVVGTRYKCSVCPDY
DLCSVCEGKGLHRGHTKLAFFSPFGHLSEGF5HSRWL
RKVKHGHFGWPGWEMGPPGNWSPRPPRAGEARPGP
TAESASGPSEDP5VNFKNVGESVAAALSPLGIEVDID
VEHGGKRSRLTPVSPES5STEKSSQPS5CCSDPSKP
GGNVEGATQ5LAEQMRKIALESEGRPEEQMESDNCS
GGDDWTHLSSKEVDPSTGELQ5LQMP5EGP5SLD
PSQEGPTGLKEAALYPHLPPEADPRLIESLSQMLSMGF
SDEGGWLRLLQTKNYDIGAALDTIQYSKHPPPL
(1-440 aa encoded by BC017222)

Reconstitution and Storage

Reconstitution:

Reconstitute at 0.25 µg/µl in 200 µl sterile water for short-term storage.

After reconstitution with sterile water, if glycerol has no effect on subsequent experiments, it is recommended to add an equal volume of glycerol for long-term storage (see Stability and Storage for more details).

If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).

Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

Shipping:

The product is shipped at ambient temperature.

Upon receipt, store it immediately at the recommended temperature (see below).

Stability and Storage

Store for up to 12 months at -20°C to -80°C as lyophilized powder.

Storage of Reconstituted Protein

Short Term Storage:

Store at 2-8°C for (1-2 weeks).

Long Term Storage:

Aliquot and store at -20°C to -80°C for up to 3 months, reconstitution with sterile water and addition of an equal volume of glycerol. Avoid repeat freeze-thaw cycles.

Selected Validation Data



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.