

**HumanKine<sup>®</sup> BMP-4 (Recombinant Human)**



Animal Component-Free	Human cell expressed	Tag-Free	Endotoxin Free
-----------------------	----------------------	----------	----------------

**Product Description**


Animal-free Recombinant Human BMP-4 is expressed in human 293 cells as a disulfide linked homodimeric glycoprotein with an apparent molecular mass of 34 kDa. This cytokine is produced in a serum-free, chemically defined media.

Alternative Names	BMP 2B, BMP 4, BMP2B, BMP2B1, BMP4, BMP-4, Bone morphogenetic protein 2B, bone morphogenetic protein 4, DVR4, MCOPS6, ZYME
Source	Human Embryonic Kidney cells (HEK293). HEK293-derived BMP-4 protein
Species Reactivity	human,mouse,bovine,Ciona intestinalis,rat


**Specifications**

Test	Method	Specification
Activity	Dose-dependent induction of alkaline phosphatase production in the ATDC-5 cell line (mouse chondrogenic cell line)	1.5-9 ng/mL EC50
Molecular Mass	SDS-PAGE	20 to 24 kDa reduced, 32 to 38 kDa non-reduced, homodimer, glycosylated
Purity	SDS-PAGE	>95%
Endotoxin	LAL	<1 EU/μg

**Activity Data**



**SDS-PAGE**



Preparation	
Shipping Temperature	ambient temperature
Formulation	2x PBS + 6% Ethanol, See Certificate of Analysis for details
Reconstitution	Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein in sterile 4 mM HCl containing 0.1% endotoxin-free recombinant human serum albumin (HSA).

Stability and Storage	Product Form	Temperature Conditions	Storage Time (From Date of Receipt)
	Lyophilized	-20°C to -80°C	Until Expiry Date
	Lyophilized	Room Temperature	2 weeks
	Reconstituted as per CofA	-20°C to -80°C	6 months
	Reconstituted as per CofA	4°C	1 week
Avoid repeated freeze-thaw cycles.			

www.ptglab.com

Document #: FR-QA118-101 Rev 0  
Data Sheet Version #: 1

**Proteintech Group, Inc.**  
5500 Pearl Street, Suite 400 Rosemont, IL 60612  
t: 1-888-478-4522; f: 1-312-455-8408  
Email: proteintech@ptglab.com