

IHC*easy* DBN1/Drebrin Ready-To-Use IHC Kit

Catalog Number: **KHC0848**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Mouse Monoclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

Drebrin (DBN1) is an actin binding and stabilizing protein with roles in endocytosis, formation of dendrite spines in neurons and coordinating cell-cell synapses in immune cells. Preferentially expressed in the brain, it is highly localized in dendritic spines and regulates spine shapes. It has two isoforms that are splice variants, with drebrin A being expressed in neural cells, and drebrin E is present in several other cell types. Recently it has been reported that circulating DBN1 levels correlated with retinal nerve fiber layer defect and may reflect the severity of retinal ganglion cells injury in glaucoma patients.


Synonyms

DoS117E, DBN1, DKFZp434D064, Drebrin, drebrin 1


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).




Immunohistochemical analysis of paraffin-embedded human gliomas tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse cerebellum tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue slide using KHC0848 (DBN1/Drebrin IHC Kit).