

# IHC*easy* SARS Ready-To-Use IHC Kit

Catalog Number: **KHC0952**

## General Information

**Sample type:**  
FFPE tissue

**Cited sample type:**

**Reactivity:**  
Human, Mouse, Rat

**Cited Reactivity:**

**Assay type:**  
Immunohistochemistry

**Primary antibody type:**  
Rabbit Polyclonal

**Secondary antibody type:**  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

## Background

SARS(Serine--tRNA ligase, cytoplasmic) Catalyzes the attachment of serine to tRNA(Ser). Is also probably able to aminoacylate tRNA(Sec) with serine, to form the misacylated tRNA L-seryl-tRNA(Sec), which will be further converted into selenocysteinyl-tRNA(Sec).


## Synonyms

SARS, Serine tRNA ligase, SERRS, SERS, seryl tRNA synthetase, Seryl tRNA(Ser/Sec) synthetase


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using KHC0952 (SARS IHC Kit).




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0952 (SARS IHC Kit).




Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0952 (SARS IHC Kit).




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0952 (SARS IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC0952 (SARS IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC0952 (SARS IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using KHC0952 (SARS IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC0952 (SARS IHC Kit).