

IHC*easy* FAM50A Ready-To-Use IHC Kit

Catalog Number: **KHC1378**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

FAM50A belongs to the FAM50 family. It is highly conserved in length and sequence across different species. It is a basic protein containing a nuclear localization signal, and may function as a DNA-binding protein or a transcriptional factor.


Synonyms

9F, DXS9928E, FAM50A, HXC 26, HXC26, Protein FAM50A, Protein HXC 26, Protein XAP 5, XAP5


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC1378 (FAM50A IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC1378 (FAM50A IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC1378 (FAM50A IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC1378 (FAM50A IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using KHC1378 (FAM50A IHC Kit).