

IHC*easy* PEG10 Ready-To-Use IHC Kit

Catalog Number: **KHC1806**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Rabbit Polyclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Synonyms

Edr, HB 1, KIAA1051, Mar2, Mart2, MEF3 like protein 1, MEF3L, MEF3L1, paternally expressed 10, PEG10, RGAG3, Ty3/Gypsy like protein


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC1806 (PEG10 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse adrenal gland tissue slide using KHC1806 (PEG10 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse skin tissue slide using KHC1806 (PEG10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC1806 (PEG10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat skin tissue slide using KHC1806 (PEG10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC1806 (PEG10 IHC Kit).