

For Research Use Only

Flow Cytometry Perm Buffer (10X)

Catalog Number: PF00011-C



www.ptglab.com

Description

Flow Cytometry Perm Buffer (10X) Provided as a concentrate which, when diluted with distilled water to a 1X solution, provide best results during intracellular staining of cytokines and other cytoplasmic antigens, by maintaining membrane permeabilization throughout staining and washing steps.

Package

150 mL


Storage

4°C

Expiration Date

6 months after receiving.

Validation Data



1X10⁶ mouse splenocytes were surface stained with CoraLite488-conjugated Anti-Mouse CD4 (GK1.5) (CL488-65104, Clone: GK1.5) and then fixed with 1X Transcription Factor Fix/Perm buffer (PF00011-A) and permeabilized with 1X Flow Cytometry Perm Buffer (PF00011-C). Cells were then stained with CoraLite647-conjugated mouse... IgG1 isotype control or 5 ul CoraLite647-conjugated Anti-Mouse Foxp3 (CL647-65089, Clone: 3G3).

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.