

For Research Use Only

Foxp3 / Transcription Factor Staining Buffer Kit

Catalog Number: PF00011



Description

The Foxp3 / Transcription Factor Staining Buffer Kit contains specially formulated buffers and solutions for optimal resolution and low background in your analysis of nuclear antigens by flow cytometry. Included are the following components, to be used during staining protocols for detection of nuclear antigens such as Foxp3 and ROR gamma.

Product Information

Components	Size
Foxp3 / Transcription Factor Fix/Perm Concentrate (4X) PF00011-A	50 mL
Foxp3 / Transcription Factor Fix/Perm Diluent (1X) PF00011-B	150 mL
Flow Cytometry Perm Buffer (10X) PF00011-C	150 mL

Package

1 kit


Storage

4°C

Expiration Date

6 months after receiving.

Validation Data



1X10⁶ mouse splenocytes were surface stained with CoraLite488-conjugated Anti-Mouse CD4 (GK1.5) (CL488-65104, Clone: GK1.5) and then fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). Cells were then stained with CoraLite647-conjugated mouse IgG1 isotype control or 5 ul CoraLite647-conjugated Ant... Mouse Foxp3 (CL647-65089, Clone: 3G3).

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.